




Microscopic and Staining Utensils




“Assistant” Micro Slide Dispenser

2440 “Assistant” Micro slide dispenser (“Objectomat DBGM”) accommodates 50 micro slides of the standard size 76 x 26 mm, dust-free storage. One micro slide always ready to hand. No finger prints by removing as slides are touched at the sides only (without slides).



1238 “Assistant” Drying rack of Plexiglas® with 12 numbered positions for drying stained micro slides. The slanted slides are set on exchangeable blotting paper. In the base of the rack is a suitable storing place for 2 packs of paper. Including 2 packs of blotting paper.

1238/1 Spare blotting paper, in packs of 25 sheets of approx. 200 x 45 mm



1240 “Assistant” Rack of Plexiglas® with numbered positions for 12 micro slides and a field for inscriptions (pencil notes can be rubbed out or removed by a moist cloth).