

“Assistent” Micro Slides

CE

“Assistent” Micro slides „ELKA“

approx. 76 x 26 mm (3 x 1")

thickness approx. 1 mm (standard thickness), selected glass quality, non-corrosive, precleaned, ready for use

50 pcs in a box,

5000 pcs in a carton

- 2400** half-white, cut edges
- 2401** half-white, ground edges, 90°
- 2402** white, cut edges
- 2403** white, ground edges, 90°

- 2404** special thickness 0,6–0,8 mm
half-white, cut edges

“Assistent” Micro slides approx. 76 x 26 mm, thickness approx. 1 mm, half-white,

with frosted end (on both sides)

- 2406** cut edges
- 2407** ground edges, 90°

with frosted end (on one side)


- 2406/1** cut edges
- 2407/1** ground edges, 90°

.../9 with tropical protection packing

Special sizes and types of micro slides are also available.



2400



2406 2407
2406/1 2407/1

“Assistent” slides Superfrost®

CE

Easy to mark and easy to read colour-coding system for medical and scientific laboratory.

The micro slides Superfrost® with colour-coded end are available in 6 different colours: white, blue, pink, yellow, green, orange. The coating is resistant to all common solvents.

- 2408** Micro slides Superfrost®, cut edges
- 2409** ground edges
- 2409/1** Micro slides, Superfrost/Plus, ground edges, with coloured label white, pack of 72 pcs.

CE



2408